Original Article

A study on the wear of enamel caused by monolithic zirconia and the subsequent phase transformation compared to two other ceramic systems

Amreen Rupawala, Smita I. Musani, Pallavi Madanshetty, Ramandeep Dugal, Umang D. Shah, Ektaa J. Sheth
Department of Prosthodontics, M.A. Rangoonwala Dental College and Research Centre, Pune, Maharashtra, India

Abstract

Purpose: The purpose of this *in vitro* research project was to evaluate and compare the wear behavior of human tooth enamel opposing monolithic zirconia and other different ceramic restorative materials and also to observe the tetragonal to monoclinic phase transformation in zirconia-based ceramics that may occur while simulating wear occurring at room temperature in a wet environment.

Materials and Methods: A total of sixty samples were prepared for this study. Fifteen discs of glazed zirconia, 15 discs of polished zirconia without glaze, 15 discs of metal ceramic, and 15 discs of lithium disilicate were fabricated. Sixty extracted premolars were collected and randomly divided into four groups of 15 each. The discs and extracted human premolars were placed onto holders on a two-body wear machine under a constant load of 5 kg to simulate the oral wear cycle. A diffractometer was used to analyze phase transformation. One-way analysis of variance and Tukey's *post hoc* tests was used.

Results: The mean loss of height of tooth samples and its standard deviation for Group I (monolithic zirconia with glaze), Group II (mechanically polished monolithic zirconia without glaze), Group III (porcelain fused to metal), and Group IV (glazed monolithic lithium disilicate) was obtained as 0.2716 ± 0.1409 , 0.1240 ± 0.0625 , 0.1567 ± 0.0996 , and 0.2377 ± 0.1350 , respectively. The highest mean loss in height was observed in Group I and the least was observed in Group II.

Conclusion: Mechanically polished zirconia showed the least amount of enamel wear followed by porcelain fused to metal and glazed monolithic lithium disilicate, whereas glazed monolithic zirconia showed the highest enamel wear.

Key Words: Lithium disilicate, phase transformation, wear, zirconia

Address for correspondence:

Dr. Ektaa J. Sheth, 702 Urvashi, Besant Street, Santacruz West, Mumbai - 400 054, Maharashtra, India. E-mail: ektaasheth@hotmail.com Received: 17th December, 2015, Accepted: 03rd May, 2016

INTRODUCTION

Ceramics are key materials for dental restorations due to their high wear resistance, biocompatibility, and above all,

Access this article online					
Quick Response Code:	Website:				
	www.j-ips.org				
100 (500) ■0000000	DOI: 10.4103/0972-4052.194940				

esthetics.^[1] Although dental ceramics have excellent properties that meet requirements of a prosthetic material, it has one major

This is an open access article distributed under the terms of the Creative Commons Attribution-NonCommercial-ShareAlike 3.0 License, which allows others to remix, tweak, and build upon the work non-commercially, as long as the author is credited and the new creations are licensed under the identical terms.

For reprints contact: reprints@medknow.com

How to cite this article: Rupawala A, Musani SI, Madanshetty P, Dugal R, Shah UD, Sheth EJ. A study on the wear of enamel caused by monolithic zirconia and the subsequent phase transformation compared to two other ceramic systems. J Indian Prosthodont Soc 2017;17:8-14.

drawback: Irreversible wear of opposing tooth structure under certain conditions. Several investigators have demonstrated that, in general, ceramic substrates cause greater and hence destructive abrasive wear of human enamel.^[2,3]

A ceramic restorative material that combines good strength without the disadvantage of increased enamel wear would be a significant addition to clinical dental practice. Although numerous in vitro wear studies have been conducted on dental ceramics, few studies have provided detailed information to characterize and predict the wear behavior of a range of dental ceramics used for all ceramic systems at present. These studies have reported the effect of individual ceramic systems on the wear of enamel.^[4,5] Taking these studies into consideration, it is challenging to compare the effect of different ceramic systems on the wear of enamel because these studies have been conducted in different environmental and testing conditions as well as in the presence of altered variables. Therefore, this in vitro research project was undertaken to evaluate and compare the wear behavior of human tooth enamel opposing different ceramic restorative materials; monolith zirconia with glaze, mechanically polished monolith zirconia without glaze, porcelain fused to metal and lithium disilicate glass ceramic using a reciprocating sliding wear testing apparatus in the same controlled environment and testing conditions using same variables.

A second part of the study is to observe the phase transformation of zirconia structure due to wear in the presence of artificial saliva. Zirconia exists in three major phases: Monoclinic, tetragonal, and cubic. The tetragonal crystalline state is responsible for the high strength and fracture toughness of yttrium oxide stabilized tetragonal zirconium dioxide polycrystals (Y-TZP) materials. Once a crack propagates within a Y-TZP material, the energy supplied by the crack can trigger the tetragonal to monoclinic phase transformation in the surrounding grains. This phase transformation leads to a local compressive stress field that hinders further crack propagation. [6,7] This so-called transformation toughening is what makes Y-TZP the strongest ceramic material class used in the dentistry. However, it is well known that on Y-TZP material surfaces, the tetragonal to monoclinic phase transformation can slowly occur by contact with water. [8-10] This is called low-temperature degradation (LTD). Zirconia has the propensity to undergo LTD, which is a kinetic phenomenon in which the polycrystalline tetragonal material slowly transforms into monoclinic zirconia over a rather narrow but important temperature range, typically between room temperature and approximately 400°C in the presence of water, depending on the stabilizer used, its concentration, and the grain size (GS) of the ceramic. This time dependant phase transformation is followed by microcracking and loss of strength which starts at the surface and then proceeds into the bulk.[11-13]

Therefore, a strict control of the phase transformation is thus imperative to retain those properties which make the ceramic material so desirable. Hence, a part of the study included observing the tetragonal to monoclinic phase transformation that may occur while simulating wear occurring at room temperature in a wet environment.

MATERIALS AND METHODS

Fabrication of zirconia discs [Box 1]

A disc of dimension 10 mm in diameter and 2 mm in height was fabricated in wax using the precision milling. Zirconia blanks were fed into the milling machine and thirty discs of the required dimensions were obtained. The discs were finished using a rubber point and later sintered. The thirty discs were randomly divided into two groups of 15 each. The discs in the first group were dried with tissue paper and glaze liquid was applied and fired. The discs of the second group were cleaned in an ultrasonic bath for 10 min and polished using Zircon-Brite polishing paste.

Fabrication of metal ceramic discs [Box 1]

A steel template with five standardized circular molds, each of diameter 10 ± 0.2 mm and depth 0.5 ± 0.1 mm was precision milled. Using this template, 15 wax patterns of standard dimensions for metal discs were fabricated and sprued. The wax pattern fabricated discs was carefully removed from the mold using a U-shaped stapler pin^[14] following the same technique used to retrieve Class II inlay wax patterns. They were invested in phosphate-bonded investment material and standardized casting procedure was carried out. The metal discs were then divested and sandblasted. The sandblasted metal discs were then secured within the molds of the refractory template and degassing was carried out. A steel template with five molds of diameter 10 ± 0.2 mm and depth 2 ± 0.1 mm was precision-milled to accommodate veneering ceramic. All the discs were layered according to manufacturer's instruction with opaque, dentin, and enamel porcelain to the desired thickness of 1.5 mm. The samples were finished and polished.

Fabrication of lithium disilicate discs [Box 1]

A steel template with five standardized circular molds, each of the diameter 10 mm and depth 2 mm, was precision-milled. Using this template, 15 wax patterns of standard dimensions were fabricated, sprued and the standardized ceramic pressing procedure was carried out. The discs were then divested. After fine divestment, the reaction layer formed during the press procedure was removed using IPS e.max Press Invex Liquid. All 15 samples were first finished using diamond burs. Zircon-Brite polishing paste was used to polish each specimen. All discs were then cleaned in an ultrasonic bath for 10 min. Discs were dried with tissue paper and glaze liquid was applied evenly and fired according to manufacturer's instructions.

Box 1: Material description

Material	Description
Zirconia discs	3M LAVA Frame Zirconia Mill Blanks (Lot number 517436) LAVA Ceram Glaze (3M LAVA, Lot number 8542C)
Metal ceramic discs	White ceramic alloy (4 All, Ivoclar Vivadent, ref. number 578942AN)
	Opaquer paste and liquid for IPS Classic (Ivoclar Vivadent, Germany, Lot P83167 and 55655)
	Feldspathic porcelain (IPS d.SIGN, Ivoclar Vivadent, Germany, dentin A1 and enamel TS-2, Lot number K34233, K3329)
Lithium disilicate discs	Lithium disilicate ingots (IPS e.Max Press, HT A1, Lot
	e.Max Ćeram Glaze (IPS e.Max, Type 1 Class 8, Lot number R81992)

All ceramic samples were then thermocycled 500 times between 5°C and 55°C.

Modified Fusayama artificial saliva was made and was stored at 37 ± 1 °C.

Sixty freshly extracted human unrestored, caries-free, nonattrited maxillary first and second premolars of young adolescent patients undergoing orthodontic extractions were collected. They were disinfected in formalin and debrided of calculus using an ultrasonic scaler and preserved in saline. The sixty extracted premolars were randomly divided into four groups of 15 each. The teeth were then mounted in autopolymerizing acrylic resin. The teeth samples were then randomly divided into four groups of 15 each:

- Group I: To abrade against glazed zirconia [Figure I]
- Group II: To abrade against polished zirconia without glaze [Figure I]
- Group III: To abrade against porcelain fused to metal [Figure I]
- Group IV: To abrade against lithium disilicate [Figure 1].

The teeth samples were placed on the worktable of the profilometry machine and the X, Y, and Z axes were adjusted. The profile of each sample of the tooth was traced and transferred to tracing paper. A perpendicular was dropped from the height of the cusp tip to the base of the tooth where it is embedded into the acrylic. This height was measured as the baseline height of that particular tooth.

The discs and extracted human premolars were placed onto holders on a two-body wear machine in which tooth sample was attached to the upper member, and the disc was attached to the lower member (rotating wheel). The disc was additionally secured with the help of condensation silicone putty. The cusp tips and ceramic discs were positioned under a constant load of 5 kg. Artificial saliva was sprayed between the tooth sample and ceramic surface at intervals during testing so as to further simulate the oral conditions. The specimens were made to rub



Figure 1: The ceramic discs Groups I, II, III, and IV

against one another in a rotating motion to simulate the oral wear cycle. The test was run for a total of 10,000 cycles on the wear machine, for each sample.

The teeth samples were then placed onto the work table of the profilometer in the same orientation as that of the first measurement. The three axes were adjusted and the height was measured again up to the base of the tooth, where it is embedded into the acrylic. A second tracing of the reduced profile of the tooth after wear was made on the previous tracing paper which contained the profile of that tooth before wear and the reduction in height was calculated.

The zirconia sample, before wear, was powdered to a fine grain and tightly packed onto the random mount present on the stage of the machine. The diffractometer was used. The knife edge attached above was first zeroed and then raised by 1.200 mm so as to contact the sample and the X-ray diffraction (XRD) was carried out. The data so obtained were assessed and plotted as a graph using OriginPro 8 (OriginLab Corporation, Northampton, Massachusetts). The graph was then compared with those in the database-PCPDF. The same protocol was carried out for the sample obtained after wear.

RESULTS

In all study groups, the height of tooth significantly reduced after the 10,000 cycles of enamel wear when compared to its initial height.

The mean loss of height of tooth samples and its standard deviation for Group I (monolithic zirconia with glaze), Group II (mechanically polished monolithic zirconia without glaze), Group III (porcelain fused to metal), and Group IV (lithium disilicate) was obtained as 0.2716 \pm 0.1409, 0.1240 \pm 0.0625, 0.1567 \pm 0.0996, and 0.2377 \pm 0.1350, respectively. The highest mean loss in height was observed in Group I and the least was observed in Group II.

Table I shows the comparison of the mean loss of height in the study groups using one-way analysis of variance. There was statistically significant difference in loss of tooth height among the study groups.

Table 2 shows the pair-wise comparison of mean loss of height in the study groups using Tukey's post hoc test. The mean difference between Group I and Group II was observed as 0.1476 with $P\!=\!0.004$. The mean difference between Group I and Group III was observed as 0.1149 with $P\!=\!0.038$. Both these are statistically significant. The mean difference between Group I and Group IV was observed as 0.0339 and was not statistically significant.

The mean difference between Group II and Group III was not statistically significant at -0.0327. The mean difference between Group II and Group IV was, however, statistically significant at -0.1137.

The mean difference between Group III and Group IV was -0.0811, with P = 0.220, hence not statistically significant.

Graph I compares the loss of height of all four groups. It is observed that Group II shows the least amount of wear while Group I shows the significantly higher wear of opposing tooth. However, Group IV shows a greater loss of height than Groups II and III but still lesser than Group I.

Zirconia specimens were analyzed on an X-ray diffractometer to observe the presence of tetragonal or monoclinic phase before as well as after wear.

Table 1: Comparison of the mean loss of height in the study groups using one-way ANOVA test

Study groups	n	Mean	SD	One-way ANOVA	
				F	P
Group I (monolithic zirconia with glaze)	15	0.2716	0.1409	5.464	0.002
Group II (monolithic zirconia without glaze but mechanically polished)	15	0.1240	0.0625		
Group III (porcelain fused to metal)	15	0.1567	0.0996		
Group IV (lithium disilicate)	15	0.2377	0.1350		


SD: Standard deviation, ANOVA: Analysis of variance

Table 2: Pair-wise comparison of mean loss of height in the study groups using Tukey's post hoc test


Group (I)	Group (J)	Mean	Significant	95% CI		
		difference (I-J)		Lower bound	Upper bound	
Group I	Group II	0.1476*	0.004	0.03746	0.25774	
	Group III	0.1149*	0.038	0.00479	0.22507	
	Group IV	0.0339	0.848	-0.07627	0.14401	
Group II	Group III	-0.0327	0.861	-0.14281	0.07747	
	Group IV	-0.1137*	0.040	-0.22387	-0.00359	
Group III	Group IV	-0.0811	0.220	-0.19121	0.02907	

^{*}The mean difference is significant at the 0.05 level. CI: Confidence interval


Graph 2 shows the XRD spectrum for zirconia specimen before wear, predominantly tetragonal phase was found in the sample.



Graph 1: Mean loss of height in all the study groups



Graph 2: X-ray diffraction spectrum of zirconia ceramic specimen from Group I before wear



Graph 3: X-ray diffraction spectrum of zirconia ceramic specimen from Group Lafter wear

Graph 3 shows the XRD spectrum for zirconia specimen after wear; this sample too showed predominantly tetragonal phase of zirconia; hence, no significant change in the phase from tetragonal to monoclinic.

DISCUSSION

Physiological wear is surface degradation that results in progressive, but the very slow loss of convexity of the cusps, which manifests as; a flattening of cusp tips on the posterior teeth and incisal edges of mammelons on the anterior teeth. However, the rate of wear may be disturbed by the introduction of restorations with wear properties differing from those of the replaced tooth structure. The severity of this problem may be accentuated by an improper occlusal scheme; for example, group function with a porcelain occlusion can cause more wear than canine-guided mutually protected occlusion. [2] Therefore, there is a need to investigate the wear properties of restorative materials and their potential abrasive effects on the opposing natural teeth.

In this *in vitro* study, wear of enamel when opposed to different ceramic systems was simulated to the oral environment using a wear machine. Two-body wear that was brought about by this machine provided a combined action of impact, followed by sliding that matches the inherent action of closure during mastication of the mandibular teeth onto the maxillary teeth for a total of 10,000 cycles.^[15] A weight of 5 kg (49 N), which is comparable to normal chewing force, was exerted onto the specimens.^[16]

Enamel varies in its properties depending on its position in the tooth and its histological structure. [17] Cuspal enamel is stronger and can withstand forces in a direction parallel to the enamel rods than perpendicular to the rods. Hence, in this study, freshly extracted nonattrited, noncarious premolars of young adolescent patients undergoing orthodontic extractions were used. Since this study was to observe the wear behavior of enamel in a clinical environment, only the cuspal tips of the dental specimens were held in contact with the ceramic specimens.

Li and Zhou^[18] investigated the influence lubrication on the wear behavior of human enamel, using a reciprocating wear test apparatus. It was found that the depth and severity of the wear scars were much smaller with artificial saliva lubrication than in dry conditions and therefore, concluded that saliva plays an important lubricant effect during the wear process of enamel in the oral environment. Because the surface hardness of glass or ceramic decreases in an aqueous environment, the two mating materials can easily adhere at the microscopic level of sharp asperities in the presence

of saliva.^[19] Consequently, a more intimate approximation between the contacting surfaces results in a higher coefficient of friction associated with exposure to aqueous media.^[20] Hence, artificial saliva was used in this study to simulate the oral environment.

When two materials in contact slide over each other, one or both of the materials will suffer wear on the surface, generally the softer of the two, i.e., tooth enamel. [21] According to this hypothesis, it was expected that being hardest of the two other ceramics under consideration, glazed zirconia as well as polished unglazed zirconia would produce the highest wear of the enamel. Our results corroborated with this hypothesis that showed that Group I (glazed monolithic zirconia) caused statistically significant wear of the opposing enamel as in the study by Heintze *et al.* where it was reported that flat glazed surfaces show more antagonistic wear than polished surfaces [22] but Group II (mechanically polished zirconia without glaze), however, showed least amount of wear of the opposing tooth enamel which was a point to ponder.

Recent studies show that a glazed zirconia surface brought about more wear than its polished surface. [23] Possible explanation is that the glazed surface is quickly worn away to reveal the rough surface of zirconia ceramic beneath. This may occur due to chairside adjustments that are made before cementation of the prosthesis or may also occur within a short period of function. The surfaces of all materials are rough at a microscopic level with sharp, rugged projections called asperities, which have a surface profile of peaks and valleys. [21] The zirconia underneath the glaze has high asperities, and owing to its high hardness value (1378–1354 Hv), tends to abrade the comparatively softer enamel opposite to it. [23]

Mechanically polished monolith zirconia (Group II) in our study showed the least amount of tooth wear among all four groups, even though it was the hardest material of the lot. Most ceramics have comparatively higher hardness values than human enamel, and the hardness of a ceramic has been used as a predictor of its potential to abrade opposing teeth.^[I] However, scientific studies have demonstrated a poor correlation between material's hardness and the abrasive potential of ceramic materials on human enamel.^[24,25] Hence, polished zirconia demonstrated less enamel wear than glazed zirconia, a finding that is in agreement with studies conducted by Preis *et al.*^[26] and Sabrah *et al.*^[27]

It was originally thought that glazed ceramics would afford the enamel a certain degree of protection against wear. However, in this study, glazed lithium disilicate (Group IV) as well as glazed porcelain fused to metal (Group III) caused greater wear of the antagonistic tooth enamel than polished unglazed zirconia.

In a study conducted by Wang et al., [28] scanning electron microscopy images of the worn surface of enamel after wear showed particles sticking on the surface of lithium disilicate antagonist as ascribed to the accumulation of chipped off enamel and glass-ceramic particles. As the lithium disilicate specimen slides over the enamel surface, sliding contact generates a frictional force which results in tensile, compressive, and shear stresses on the enamel surface. Microcracks develop within the subsurface and these subsurface cracks propagate to eventually form a particle which becomes displaced according to Reid et al. [29]

As stated by Seghi et al., [24] Dicor glass ceramic was found to be significantly harder than other glass ceramics but had lower relative abrasiveness. Lack of significant crystalline phase may account for this lower abrasiveness despite its high hardness value. However, pressable lithium disilicate material consists of higher crystalline content (70% volume), which may be a contributing factor for increased abrasiveness of antagonistic tooth enamel, as seen in this study.

As evidenced by the literature and the results of this study, it may be prudent to conclude that polishing the surface of the ceramic is the wisest choice of surface finish when opposing natural dentition.

Powder XRD analysis was done on the zirconia specimens before and after wear simulation on Bruker D8 Advance X-ray diffractometer using 2θ – θ method, where θ is the angle of reflection. The samples were observed for tetragonal and monoclinic phases of zirconia and the data so collected were interpreted to achieve the zirconia spectra graph. In the spectra obtained by the 2θ – θ method in the first zirconia sample before wear, strong peaks representing the tetragonal phase were observed, while peaks assigned to the monoclinic phase were not distinct. In the XRD of the sample after wear in the presence of artificial saliva, peaks of tetragonal phase were observed only, with no presence of transformation-produced monoclinic phase. Irrespective of the mechanism, it is evident that the presence of moisture as well as increased temperatures can induce this t \rightarrow m transformation.

The XRD results of the present study are in agreement with the results of the study carried out by Eástková et al. [30] in 2004 in which absolutely no layer of transformed zirconia was recorded when specimen tested at $140^{\circ}\text{C}/15$ h. The transformed layer was found at a temperature of 180°C and higher only, the extent of $t \to m$ transformation increasing with the rise in temperature. Numerous researchers have reported that reducing the average GS in zirconia-based ceramics has a beneficial effect on the stability of the tetragonal phase [31,32] and therefore on LTD. The absence of transformation of phase in the specimen used in this study may be due to the reduced GS of the material,

 $0.5 \, \mu m$ which is well within the limits of standard GS allowed for zirconia intended for dental use.

Thus, taking all these factors into consideration, the operator must be capable of making an informed decision when selecting a restorative material which will restore form, function, and esthetics as well as protect the present dentition from additional damage.

CONCLUSION

- Mechanically polished zirconia showed the least amount
 of wear of enamel as compared to glazed monolithic
 zirconia thus suggesting that polishing the monolithic
 zirconia surface might be the best treatment to reduce
 surface roughness and antagonistic wear of enamel. This
 is beneficial in high-load bearing areas, making polish a
 superior choice over an over-glaze
- Lithium disilicate showed statistically significant amount of antagonistic enamel wear thus making it a poor choice of prosthesis in patients with parafunctional habits
- Glazed porcelain fused to metal caused greater wear of the antagonistic tooth enamel than polished unglazed zirconia
- The zirconia showed almost no t → m transformation after 10,000 cycles of wear in a wet environment thus preventing any detrimental changes in its mechanical properties when used in the oral cavity.

Financial support and sponsorship Nil.

Conflicts of interest

There are no conflicts of interest.

REFERENCES

- Kelly J. Ceramics in restorative and prosthetic dentistry. Annu Rev Mater Sci 1997;27:443-68.
- Wiley MG. Effects of porcelain on occluding surfaces of restored teeth. J Prosthet Dent 1989;61:133-7.
- 3. Mahalick JA, Knap FJ, Weiter EJ. Occusal wear in prosthodontics. J Am Dent Assoc 1971;82:154-9.
- Mulay G, Dugal R, Buhranpurwala M. An evaluation of wear of human enamel opposed by ceramics of different surface finishes. J Indian Prosthodont Soc 2015;15:111-8.
- Khandelwal M, Jain D. A comparative evaluation of wear of enamel caused by ceramics with different fusion temperatures. J Indian Prosthodont Soc 2013;13:513-9.
- Kelly JR, Denry I. Stabilized zirconia as a structural ceramic: An overview. Dent Mater 2008;24:289-98.
- Lughi V, Sergo V. Low temperature degradation aging-of zirconia: A critical review of the relevant aspects in dentistry. Dent Mater 2010;26:807-20.
- Chevalier J, Gremillard L, Deville S. Low-temperature degradation of zirconia and implications for biomedical implants. Annu Rev Mater Res 2007;37:1-32.
- Chevalier J, Gremillard L, Virkar AV, Clarke DR. The tetragonal-monoclinic transformation in zirconia: Lessons learned and future trends. J Am Ceram Soc 2009;92:1901-20.

- Chevalier J. What future for zirconia as a biomaterial? Biomaterials 2006:27:535-43.
- Ban S, Suehiro Y, Nakanishi H, Nawa H. Fracture toughness of dental zirconia before and after autoclaving. J Ceram Soc Jpn 2010;118:406-9.
- Eichler J, Rodel J, Ulrich E, Mark H. Effect of grain size on mechanical properties of submicrometer 3Y-TZP: Fracture strength and hydrothermal degradation. J Am Ceram Soc 2007;90:2830-6.
- Lin JD, Duh JG. Crystallite size and microstrain of thermally aged low-ceria-and low-yttria-doped zirconia. J Am Ceram Soc 1998;81:853-60.
- Soratur SH. Essentials of Prosthodontics. New Delhi: Jaypee Brothers Publishers; 2009. p. 159.
- Elmaria A, Goldstein G, Vijayaraghavan T, Legeros RZ, Hittelman EL. An evaluation of wear when enamel is opposed by various ceramic materials and gold. J Prosthet Dent 2006;96:345-53.
- Jung YS, Lee JW, Choi YJ, Ahn JS, Shin SW, Huh JB. A study on the *in-vitro* wear of the natural tooth structure by opposing zirconia or dental porcelain. J Adv Prosthodont 2010;2:111-5.
- Phillips RW. Skinners Science of Dental Materials. 8th ed. Philadelphia: W.B Saunders Co., Harcourt Brace and Co.; 1982.
- Li H, Zhou Z. Wear behaviour of human teeth in dry and artificial saliva conditions. Wear 2002;249:980-4.
- Callister WD Jr. Materials Science and Engineering: An Introduction. 4th ed. New York: John Wiley & Sons; 1997. p. 179-95.
- Koran A, Craig RG, Tillitson EW. Coefficient of friction of prosthetic tooth materials. J Prosthet Dent 1972;27:269-74.
- Hutchings I. Tribology: Friction and Wear of Engineering Materials. London: Butterworth-Heinemann Ltd.; 1992.
- 22. Heintze SD, Cavalleri A, Forjanic M, Zellweger G, Rousson V. Wear of

- ceramic and antagonist A systematic evaluation of influencing factors in vitro. Dent Mater 2008;24:433-49.
- Janyavula S, Lawson N, Cakir D, Beck P, Ramp LC, Burgess JO. The wear of polished and glazed zirconia against enamel. J Prosthet Dent 2013;109:22-9.
- Seghi RR, Rosenstiel SF, Bauer P. Abrasion of human enamel by different dental ceramics in vitro. J Dent Res 1991;70:221-5.
- Won-Suck O, Delong R, Anusavice KJ. Effect of restorative materials on wear of human enamel. J Prosthet Dent 2002;64:194-203.
- Preis V, Behr M, Kolbeck C, Hahnel S, Handel G, Rosentritt M. Wear performance of substructure ceramics and veneering porcelains. Dent Mater 2011;27:796-804.
- Sabrah AH, Cook NB, Luangruangrong P, Hara AT, Bottino MC. Full-contour Y-TZP ceramic surface roughness effect on synthetic hydroxyapatite wear. Dent Mater 2013;29:666-73.
- Wang L, Liu Y, Wenjie S, Feng H, Tao Y, Ma Z. Friction and wear behaviors of dental ceramics against natural tooth enamel. J Eur Ceram Soc 2012;32:2599-606.
- Reid CN, Fisher J, Jacobsen PH. Fatigue and wear of dental materials. J Dent 1990;18:209-15.
- Eástková K, Hadraba H, Cihlá J. Hydrothermal ageing of tetragonal zirconia ceramics. Ceram Silikáty 2004;48:85-92.
- Garvie RC. Stabilization of the tetragonal structure in zirconia microcrystals. J Phys Chem 1978;82:218-24.
- Tsukuma K, Kubota Y, Tsukidate T. Thermal and mechanical properties of Y2O3-stabilized tetragonal zirconia polycrystals. In: Claussen N, Ruehle M, Heuer AH, editors. Science and Technology of Zirconia II. Columbus, OH: The American Ceramic Society; 1984. p. 382-90.



"Quick Response Code" link for full text articles

The journal issue has a unique new feature for reaching to the journal's website without typing a single letter. Each article on its first page has a "Quick Response Code". Using any mobile or other hand-held device with camera and GPRS/other internet source, one can reach to the full text of that particular article on the journal's website. Start a QR-code reading software (see list of free applications from http://tinyurl.com/yzlh2tc) and point the camera to the QR-code printed in the journal. It will automatically take you to the HTML full text of that article. One can also use a desktop or laptop with web camera for similar functionality. See http://tinyurl.com/2bw7fn3 or http://tinyurl.com/3ysr3me for the free applications.